1539644.9

Instructions for Use

Cordis OPTEASE® Retrievable Vena Cava Filter

Mode d'emploi

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
English

STERILE. Sterilized with ethylene oxide gas. Nonpyrogenic. Radiopaque.
FOR SINGLE USE ONLY. DO NOT RESTERILIZE.

I. Device Description

The Cordis **OPTEASE®** Retrievable Vena Cava Filter (Retrieval Filter) is designed for percutaneous delivery of a retrievable vena cava filter to the inferior vena cava (IVC).

The self-centering **OPTEASE®** Retrievable Filter is laser cut from nickel titanium alloy (Nitinol) tubing. The proximal and distal baskets of the **OPTEASE®** Retrieval Filter, which consist of struts in a six diamond shape configuration, are designed for optimal dot capture. The baskets are connected by six straight struts. A single row of fixation barbs is present at the cranial end of the struts. These barbs, intended for fixation to the vessel wall, are extensions of the parallel struts. A hook is centrally located at the caudal basket extremity and allows for filter retrieval using a snare (Figure 1).

The constrained filter is flexible and achieves its unconstrained diameter upon deployment in the IVC. Upon deployment, the filter imparts an outward radial force on the luminal surface of the vena cava to ensure proper positioning and stability. The **OPTEASE®** Retrieval Filter is designed to prevent pulmonary embolism while maintaining caval patency.

The introduction kit is comprised of a filter storage tube, a 6F (2.0 mm ID) Cordis **BRITE TIP®** Catheter Sheath Introducer (Sheath Introducer), an Angiographic Vessel Dilator with an open end, side holes and two marker bands at a distance of 30 mm (end-to-end), and an obturator with marker and storage tube securement bump.

The Angiographic Vessel Dilator is designed to also provide angiographic visualization and linear measurement of the IVC when combined with the delivery of radiopaque contrast media.

The obturator serves to advance the filter from the storage tube into the Sheath Introducer and to advance the filter through the Sheath Introducer to the implantation site. The constrained filter is supplied in a plastic storage tube, which is to be loaded as a system into the Sheath Introducer hemostasis valve. The storage tube is printed with colored arrows and text (femoral: green; jugular/antecubital: blue) to indicate the correct orientation. The arrow of the desired access site will point into the Sheath Introducer hemostasis valve.

The length of the introduction kit, the recommended guidewire length and diameter, and the recommended venous access sites are given in Table I.

Table I.

Catalog Number	Usable BRITE TIP® Sheath Introducer Length (cm)		Angiographic Vessel Dilator Length (cm)
466-F210A	55		66
466-F210B	90		101
Catalog Number	Recommended Guidewire		Recommended Venous Access Sites
	Length (cm)	Diameter (in.)	
466-F210A	150	0.035" (0.89 mm) Cat. No. 502-521	Femoral & Jugular
466-F210B	260	0.035" (0.89 mm) Cat. No. 502-455	Femoral, Jugular & Antecubital

Retrieval of the **OPTEASE®** Retrieval Filter is from femoral venous access only (refer to Section VIII, Table III).

II. Indications

The **OPTEASE®** Retrieval Filter is indicated for the prevention of pulmonary embolism (PE) via percutaneous placement in the IVC in patients considered at high risk of PE.

The Angiographic Vessel Dilator is designed to provide angiographic visualization and linear measurement of the vasculature when combined with the delivery of radiopaque contrast media to the vena cava.

III. Contraindications

Vena cava filters should not be implanted in:

- Patients with risk of septic embolism.
- Patients with uncontrolled infectious disease.
- Patients with an IVC diameter larger than 30 mm.
- Patients contraindicated for procedures under fluoroscopy.
- Patients with demonstrated hypersensitivity to one of the components of the **OPTEASE®** Retrieval Filter.

There are no known contraindications for use of the Angiographic Vessel Dilator.

IV. Warnings

- This product is designed and intended for single use. It is not designed to undergo reprocessing and re-sterilization after initial use. Reuse of this product, including after reprocessing and/or re-sterilization, may cause a loss of structural integrity which could lead to a failure of the device to perform as intended and may lead to a loss of critical labeling/use information all of which present a potential risk to patient safety.
- The **OPTEASE®** Retrieval Filter has not been studied for long term implantation and must be retrieved within 12 days after placement.
- Retrieval of the **OPTEASE®** Retrieval Filter is possible only from femoral vein approach. Retrieval is performed using the Cordis **OPTEASE®** Retrieval Catheter (catalog number, 466-C210F) and the recommended accessories (refer to Section VIII, Table III). These devices used for filter retrieval are not included in the **OPTEASE®** Retrieval Filter Introduction Set. The IFU for **OPTEASE®** Retrieval Filter retrieval is contained in the Cordis **OPTEASE®** Retrieval Catheter packaging.
- When injecting contrast medium through the Angiographic Vessel Dilator, do not exceed maximum pressure rating of 800 psi. Ensure that a high-pressure connection line is used.
- After filter implantation, any catheterization procedure requiring passage of a device through the filter may be impeded.
- The **OPTEASE®** Retrieval Filter is supplied constrained in a plastic storage tube indicating the appropriate orientation for femoral and jugular/antecubital approaches. **Never reload a fully ejected filter into the storage tube as this could result in incorrect orientation for the selected access site.** Never reload a (partially) ejected filter into the storage tube as this could affect its shape, and function. Accordingly, Cordis will not be responsible for any direct, incidental or consequential damages resulting from replacement of the **OPTEASE®** Retrieval Filter in the plastic storage tube.
- Filter fracture is a potential complication of vena cava filters. Anatomic locations that create concentrated stress points from filter deformation (for example, deployment at apex of scoliosis, overlapping of either of the renal ostia, or placement adjacent to a vertebral osteophyte) may contribute to fracture of a particular filter strut. However, reports of adverse clinical sequelae from filter fractures are rare.
- The **OPTEASE®** Retrieval Filter should only be used by physicians who are trained in diagnostic and percutaneous interventional techniques, for instance placement of vena cava filters. Accordingly, Cordis will not be responsible for any direct or consequential damages or expenses resulting from use by untrained personnel.
- Persons with allergic reactions to nickel titanium (Nitinol) may suffer an allergic response to this implant.
- Retrieval of the **OPTEASE®** Retrieval Filter with a filter fracture present may result in complications.
- Persons with allergic reactions to nickel may suffer an allergic response to this device.

V. Precautions

- Store in a cool, dark, dry place.
- Do not use if package is open or damaged.
- Use prior to "Use By" date.
- Do not resterilize by any other method. Do not expose to temperatures above 54°C (130°F) for any length of time.
- Do not expose to organic solvents.
- The **OPTEASE®** Retrieval Filter has been tested and qualified with the accompanying accessories. The use of any other accessory could result in complications and/or an unsuccessful procedure.
- If strong resistance is met during any stage of the procedure, discontinue the procedure and determine the cause before proceeding.

MRI Safety Information

A patient with the **OPTEASE®** Vena Cava Filter may be safely scanned under the following conditions.
Failure to follow these conditions may result in injury to the patient.

Name/identification of device	Cordis OPTEASE® Vena Cava Filters
Nominal values of static magnetic field (T)	1.5 T and 3.0 T
Maximum spatial field gradient (T/m) and (Gauss/cm)	25 T/m (2500 Gauss/cm)
RF excitation	Circularly polarized (CP)
RF transmit coil type	Whole body transmit coil Head RF transmit-receive coil
RF receive coil type	Any receive only coil may be used
Maximum whole body SAR (W/kg)	4.0 W/kg
Limits on scan duration	15 minutes of continuous RF (a sequence or back to back series/scan without breaks) followed by a wait time of 10 minutes if this limit is reached
MR image artifact	The presence of this implant produced an image artifact of approximately 9 mm when imaged with a gradient echo pulse sequence and a 3.0 T MRI system
If information about a specific parameter is not included, there are no conditions associated with that parameter	

VI. Complications

Procedures requiring percutaneous interventional techniques should not be attempted by physicians unfamiliar with the possible complications. Complications may occur at any time during the procedure.

Possible procedure complications include, but are not limited to, the following:

- air embolism
- hematoma at the puncture site
- incorrect positioning of the filter
- incorrect orientation of the filter

- perforation of the vessel wall
- restriction of blood flow
- occlusion of small vessel
- distal embolization
- infection
- intimal tear
- filter fracture
- thrombus formation.

If the filter is not removed within 12 days of implantation, the possible long-term complications associated with filter implantation include, but are not limited to the following:

- filter obstruction
- filter perforation of the vena cava wall
- filter migration
- filter fracture
- recurrent pulmonary embolism.

Filter fracture is a potential complication of vena cava filters. Anatomic locations that create concentrated stress points from filter deformation (for example, deployment at apex of scoliosis, overlapping of either of the renal ostia, or placement adjacent to a vertebral osteophyte) may contribute to fracture of a particular filter strut. However, reports of adverse clinical sequelae from filter fractures are rare.

Retrieval of the Cordis **OPTEASE®** Retrievable Filter with a filter fracture present may result in complications.

VII. Recommended Percutaneous Procedure for Filter Implantation

1. Select a suitable venous access site.
Note: It may be preferable to use the right-sided approach because filter advancement can be problematic when using a tortuous left-sided vessel approach.
2. Remove the components of the introduction kit from the package using sterile technique.
3. Remove air from the Sheath Introducer by flushing with heparinized saline or suitable isotonic solution.
4. Insert the Angiographic Vessel Dilator through the Sheath Introducer hemostasis valve, snapping it into place at the hub (Figure 3). Flush with heparinized saline or suitable isotonic solution.
5. A Cordis J-tip 0.035" (0.89 mm) **EMERALD®** guidewire is recommended to introduce the Sheath Introducer and the Angiographic Vessel Dilator. The length and the maximum diameter of the recommended guidewire are listed in Table I.
6. Introduce the cannula of an angiographic needle into the vessel using aseptic technique.
7. Slide the guidewire introducer over the J tip of the guidewire to straighten it for insertion. Holding the needle in place, insert the guidewire through the needle and into the vessel. Gently advance the guidewire to the desired location.
8. Holding the guidewire in place, withdraw the needle and apply pressure to the puncture site until the Sheath Introducer with the Angiographic Vessel Dilator is inserted into the vasculature. Pass the Sheath Introducer system over the guidewire, grasping the Sheath Introducer close to the skin to prevent buckling. Using a rotating motion, advance the assembly through the tissue and into the vessel (Figure 4).
9. Position the Sheath Introducer radiopaque tip and the marker bands of the Angiographic Vessel Dilator in the inferior vena cava below the renal veins in preparation for an angiographic overview of the IVC.
10. Remove the guidewire.
11. Determine the diameter of the inferior vena cava at the intended implantation site below the lowest renal vein by injecting contrast medium through the Angiographic Vessel Dilator and using its marker bands as a reference. The distance between the two marker bands, end-to-end, is 30 mm.
Caution: Do not use with Ethiodol or Lipiodol® contrast media, or other such contrast media, which incorporate the components of these agents.
Caution: Do not exceed 800 psi when injecting through a power injector. Ensure that a high-pressure connection line is used.
12. Reintroduce the guidewire. Advance the Sheath Introducer tip to the desired location in the IVC.
13. Remove the guidewire. Detach the Angiographic Vessel Dilator from the Sheath Introducer by releasing the snap-fit ring at the hub (Figure 5). Withdraw the Angiographic Vessel Dilator.
Caution: To avoid damage to the Sheath Introducer tip, do not withdraw the dilator until the Sheath Introducer tip is at the desired location in the IVC.
Caution: Do not leave a Sheath Introducer in place for extended periods of time without a dilator or an obturator to support the cannula wall.
14. Aspirate from the sideport extension to remove any potential air.
15. According to the selected venous access site, determine which end of the storage tube (containing the filter) is to be placed into the Sheath Introducer hemostasis valve. This is indicated by the printed colored arrows and text (**femoral: green; jugular/ antecubital: blue**) on the storage tube. *The arrow of the desired access site will point into the Sheath Introducer hemostasis valve* (Figure 6A/B).
16. Place the appropriate end of the storage tube (containing the **OPTEASE®** Retrievable Filter), as far as possible into Sheath Introducer hemostasis valve (Figure 6A/B).
17. Slowly advance the filter into the Sheath Introducer by advancing the obturator through the end of the storage tube until the filter is positioned well into the cannula of the Sheath Introducer.
18. To restore hemostasis, withdraw the storage tube out of the Sheath Introducer hemostasis valve. The storage tube can be fixed on the hub-end of the obturator by sliding it over the bump (figure 7A/B).
19. Continue to advance the filter until the marker on the obturator is positioned at the Sheath Introducer hemostasis valve. This indicates that the filter is at the distal tip of the Sheath Introducer but still fully within the Sheath Introducer (Figure 7A/B).
Note: If filter advancement is problematic when using a tortuous vessel approach, stop filter advancement prior to the curve. Advance the Sheath Introducer to negotiate the curve, then continue to advance the filter.
20. Filter release (or deployment) is performed under continuous fluoroscopy (Figure 8A/B). Prior to releasing the **OPTEASE®** Retrievable Filter from the Sheath Introducer ensure that:
 - a. the intended filter location in the inferior vena cava is correct. Should the filter location be incorrect, reposition the Sheath Introducer; and
 - b. the filter retrieval hook is orientated towards the caudal end of the vena cava. Should the filter orientation be incorrect, remove the Sheath Introducer (with filter inside) from the patient and discard the system. Recommence the procedure using a new sterile Sheath Introducer and storage tube with filter.**Note:** The constrained, unexpanded, length of the **OPTEASE®** Retrievable Filter is 67 mm. During release from the Sheath Introducer, the length of the filter will shorten as it expands. The ends of the filter will shorten equally as they move toward the axial center of the filter. Filter shortening is depicted in Figure 2. The length of the deployed filter is determined by its expanded diameter in the IVC. The shortening at each end of the expanded filter and the corresponding filter length is provided in Table II. The position of the expanded filter within the IVC can be determined by selecting either the midpoint (i.e., axial center) or one end of the constrained filter as a reference point. When the filter is positioned by using the midpoint of the constrained filter as a reference point, the position of the midpoint of the expanded filter will not change. When the filter is positioned by selecting one end of the constrained filter as a reference point, the end of the expanded filter will move (i.e., filter shortening) toward the axial center.

Table II. Filter Length per IVC Diameter
Unexpanded Filter Length = 67 mm

IVC Diameter (mm)	Expanded Filter Length (mm)	Filter Shortening at each end (mm)
10	66	0.5
12	65	1
14	64	1.5
16	64	1.5
18	63	2
20	62	2.5
22	61	3
24	59	4
26	57	5
28	56	5.5
30	54	6.5

Note: Filter length and shortening values given are approximate.

21. To release the filter, keep the obturator fixed in position, then pull the Sheath Introducer back over the obturator (Figure 8A/B).
22. Ensure that the filter is fully released and deployed.
23. Once the **OPTEASE®** Retrievable Filter has been deployed, remove the obturator from the Sheath Introducer.
24. Perform a control cavogram before terminating the procedure.
25. Remove the Sheath Introducer when clinically indicated by placing compression on the vessel above the puncture site, and slowly withdrawing the Sheath Introducer.
26. Discard the introduction kit and packaging materials.
Note: After use, the introduction kit and packaging materials may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and with applicable local, state and federal laws and regulations.

VIII. Recommended Devices Required for the Retrieval of the Cordis **OPTEASE®** Retrievable Filter

The **OPTEASE®** Retrievable Filter has not been studied for long term implantation and must be retrieved within 12 days after placement.

The recommended devices required for the retrieval of the **OPTEASE®** Retrievable Filter are listed in Table III below.

The IFU for the retrieval of the **OPTEASE®** Retrievable Filter is contained in the packaging of the Cordis **OPTEASE®** Retrieval Catheter.

Table III. Recommended Devices for the Retrieval of the **OPTEASE®** Retrievable Filter

Recommended Devices	Dimension	Cat. No.
OPTEASE® Retrieval Catheter with radiopaque tip	10F, 80 cm length	466-C210F
Cordis J-tip EMERALD® 0.035" (0.89 mm) Guidewire	260 cm length	502-455
Cordis AVANTI® + Catheter Sheath Introducer	10F (3.3 mm), 11 cm length	504-610X
Goose Neck Snare Kit	10 mm loop diameter, 4F catheter 15 mm loop diameter, 6F catheter 20 mm loop diameter, 6F catheter 25 mm loop diameter, 6F catheter	GN1000 GN1500 GN2000 GN2500

* Ethiodol and Lipiodol are trademarks of Guerbet S.A.

30 mm loop diameter, 6F catheter
All 120 cm snare length and 102 cm catheter length

GN3000

IX. Disclaimer of Warranty and Limitation of Remedy
There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the Cordis product(s) described in this publication. Under no circumstances shall Cordis be liable for any direct, incidental, or consequential damages other than as expressly provided by specific law. No person has the authority to bind Cordis to any representation or warranty except as specifically set forth herein.

Descriptions or specifications in Cordis printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties.

Cordis Corporation will not be responsible for any direct, incidental, or consequential damages resulting from reuse of the product.

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© 2023 Cordis. All rights Reserved. CORDIS, the Cordis LOGO, OPTEASE, BRITE TIP, EMERALD and AVANTI are trademarks of Cordis and may be registered in the US and/or in other countries. All other marks are the property of their respective owners.
* Ethiodol i Lipiodol trgovački su znakovi tvrtke Guerbet S.A.